


# Newton's Law of Cooling and the Heat Transfer

PhyPiDAQ  
Digital Measurement System Based on  
Raspberry Pi



## Objectives:

- Measure the temperature of a cooling fluid against time by means of a 100Pt resistance temperature detector (RTD) with MAX31865 sensor connected at the Raspberry Pi.
- Use various capabilities of the PhyPiDAQ-Software to visualise the exponentially decreasing process of the temperature difference of a cooling body with respect to the surroundings in real-time.
- Employ spreadsheets like LibreOffice or Excel to calculate the decreasing rate of temperature for different substances in various quantities.



*Various graphical representations on the PhyPiDAQ-window as measuring the temperature against time in cooling processes for different substances.*

## Configurations:

-Configure the experiment and the MAX31865 sensor on the Graphical Interface of the PhyPiDAQ Software according to  
`Temperatur_Max31865_wather.daq`  
`MAX31865Config.yaml`

## Setup

