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```
# example of a configuration file for ADC ADS1115

DAQModule: ADS1115Config

ADCChannels: [0]    # active ADC-Channels
                    # possible values: 0, 1, 2, 3
                    # when using differential mode:
                    # - 0 = ADCChannel 0
                    #   minus ADCChannel 1
                    # - 1 = ADCChannel 0
                    #   minus ADCChannel 3
                    # - 2 = ADCChannel 1
                    #   minus ADCChannel 3
                    # - 3 = ADCChannel 2
                    #   minus ADCChannel 3
DifModeChan: [false] # enable differential mode for Channels
Gain: [2/3]          # programmable gain of ADC-Channel
                    # possible values for Gain:
                    # - 2/3 = +/-6.144V
                    # - 1 = +/-4.096V
                    # - 2 = +/-2.048V
                    # - 4 = +/-1.024V
                    # - 8 = +/-0.512V
                    # - 16 = +/-0.256V
sampleRate: 860     # programmable Sample Rate of ADS1115
                    # possible values for SampleRate:
                    # 8, 16, 32, 64, 128, 250, 475, 860
```

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Abb.13

```
#
# -- configuration options for Channels
#
#ChanLabels: ['X1', 'X2'] # names for channels
#ChanUnits: ['a.U.', 'a.U.'] # units for channels
#ChanLabels: [U, U] # names for channels
ChanUnits: [V] # units for channels
ChanNams: [U_G]
ChanColors: [darkblue, sienna] # channel colours in display

# eventually overwrite Channel Limits obtained from device config
ChanLimits:
- [0., 4.] # chan 0
## - [0., 1.] # chan 1
## - [0., 1.] # chan 2

# calibration of channel values
# - null or - <factor> or - [ [ <true values> ], [ <raw values> ] ]
#ChanCalib:
# - 1. # chan0: simple calibration factor
# - [ [0.,1.], [0., 1.] ] # chan1: interpolation: [true]([<raw>] )
# - null # chan2: no calibration

# apply formulae to (calibrated) channel values
#ChanFormula:
# - c0 + c1 # chan0
# - c1 # chan1
# - null # chan2 : no formula
```

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```
#
# -- configuration options for graphical display
#
Interval: 1.0          # logging interval
#NHistoryPoints: 120   # number of points used in history buffer
#DisplayModule: DataLogger # history of channel signals
DisplayModule: DataGraphs # text, bar-graph, history and xy-view
Title: "Photoeffekt546nm" # display title
#XYmode: false        # enable/disable XY-display
## if more than two channels active:
#Chan2Axes: [0, 1, 0]  # assign channels to axes
#xyPlots:              # define which axes to show
# - [0, 1]             # in xy-plot
# - [0, 2]
# - [1, 2]

#
# -- configuration options for output to file
#
DataFile: 546nm.csv    # file name for output file,
#DataFile: null       # null to disable
CSVseparator: ';'     # field separator, set to ';' for German Excel

# enable buffering of latest data (depth NHistoryPoints from above)
#bufferData: PhyPiData # file name to track latest data and eventually
#bufferData: null      # store them, or null to switch off

# control status LEDs
#RunLED: 20 # display run status on GPIO pin 20
#ReadoutLED: 21 # display readout on GPIO pin 21
```

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